

Phospho-Irs1-S307(Human Specific) Rabbit Pab For Wb, Ihc-P, Elisa - P35568

Item Number: CM0001499



Introduction

Phospho-specific rabbit polyclonal antibody targeting IRS1 pS307 (human specific), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for insulin signaling, metabolic disease research, and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Phospho-IRS1-S307(Human Specific) Rabbit pAb
Remarks/Alias	HIRS-1; Phospho-IRS1-S307(Human Specific)
Species	Human
GeneID	3667
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S307 of human IRS1 (NP_005535.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35568
Protein Weight	132kDa
Shipping	Ice bag