

Phospho-Irs1-S312 Rabbit Polyclonal Antibody For Western Blot - P35568

Item Number: CM0001213



Introduction

Phospho-specific rabbit polyclonal antibody targeting human IRS1 at serine 312. Validated for Western blotting with cross-reactivity to mouse and rat, suitable for insulin signaling and metabolism studies. - P35568

[Learn More](#)

Parameter	Value
Product Name	Phospho-IRS1-S312 Rabbit pAb
Remarks/Alias	HIRS-1; Phospho-IRS1-S312
Species	Human
Gene ID (Human)	3667
Gene ID	3667
Immunogen	Synthetic peptide - A phospho specific peptide corresponding to residues surrounding S312 of human IRS1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human, Mouse, Rat
SWISS	P35568
Protein Weight	132kDa
Shipping	Ice bag