

Phospho-Irs1-S612/S636 Rabbit Mab For Wb, Ihc-P, Elisa - P35568

Item Number: CM0000749



Introduction

Phospho-IRS1-S612/S636 rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Reacts with human and mouse IRS1. Key tool for insulin signaling and metabolic disorder studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-IRS1-S612/S636 Rabbit mAb
Remarks/Alias	HIRS-1; Phospho-IRS1-S612/S636
Species	Human
GeneID (Human)	3667
GeneID	3667
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S616 of human IRS1 (NP_005535.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P35568
Protein Weight	132kDa
Shipping	Ice bag