

Phospho-Irs1-S636/639 Rabbit Pab For Wb, Elisa - P35568

Item Number: CM0001602



Introduction

Phospho-specific rabbit polyclonal antibody targeting IRS1 phosphorylated at Ser636 and Ser639. Validated for Western blot and ELISA in human samples. Explore insulin signaling and diabetes research with this high-quality antibody.

[Learn More](#)

Product Name	Remarks/Alias	Species	GenelD (Human)	GeneID	Immunogen	Source	Category	Application	Cross Reactivity	SWISS	Protein Weight	Shipping
Phospho-IRS1-S636/639 Rabbit pAb	HIRS-1; Phospho-IRS1-S636/639	Human	3667	3667	Synthetic peptide: A synthetic phosphorylated peptide around S636 & S639 of human IRS1 (NP_005535.1).	Rabbit	Phosphorylated Antibodies	WB, ELISA	Human	P35568	132kDa	Ice bag