

Phospho-Jak2-Y1007/1008 Rabbit Mab For Wb, Elisa - O60674

Item Number: CM0001440



Introduction

Phospho-JAK2-Y1007/1008 rabbit monoclonal antibody specifically detects phosphorylated human JAK2 at Y1007/Y1008. Validated for WB and ELISA, ideal for studying JAK-STAT signaling in cancer and immunity research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-JAK2-Y1007/1008 Rabbit mAb
Remarks/Alias	JTK10
Species	Human
GeneID (Human)	3717
GeneID	3717
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y1007 & Y1008 of human JAK2 (NP_001309123.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O60674
Protein Weight	131kDa
Shipping	Ice bag