

Phospho-Jak2-Y1008 Rabbit Pab For Wb, Elisa - O60674

Item Number: CM0001758



Introduction

Phospho-JAK2-Y1008 rabbit polyclonal antibody validated for WB and ELISA. Detects human JAK2 phosphorylated at Tyr1008, a critical activation site in cytokine signaling and chromatin regulation. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-JAK2-Y1008 Rabbit pAb
Remarks/Alias	JTK10; Phospho-JAK2-Y1008
Species	Human
GeneID (Human)	3717
GeneID	3717
Immunogen	A synthetic phosphorylated peptide around Y1008 of human Phospho-JAK2-Y1008 (NP_004963.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O60674
Protein Weight	131kDa
Shipping	Ice bag