

Phospho-Jak2-Y221 Rabbit Polyclonal Antibody For Wb And Elisa - O60674

Item Number: CM0000913



Introduction

Phospho-JAK2-Y221 rabbit polyclonal antibody detects endogenous phosphorylated JAK2 at tyrosine 221 in human, mouse, and rat. Ideal for WB and ELISA applications in JAK/STAT signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-JAK2-Y221 Rabbit pAb
Remarks/Alias	JTK10; Phospho-JAK2-Y221
Species	Human
GeneID (Human)	3717
Immunogen	A synthetic phosphorylated peptide around Y221 of human JAK2 (NP_004963.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60674
Protein Weight	131kDa
Shipping	Ice bag