

Phospho-Jnk1-T183/Y185 + Jnk2-T183/Y185 + Jnk3-T221/Y223 Rabbit Pab - P45983 / P45984 / P53779

Item Number: CM0000687



Introduction

Phospho-specific rabbit polyclonal antibody targeting JNK1/2/3 activation sites, validated for Western blot and ELISA. Designed for human and mouse research in stress signaling and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-JNK1-T183/Y185 + JNK2-T183/Y185 + JNK3-T221/Y223 Rabbit pAb
Species	Human
GeneID (Human)	5599/5601/5602
GeneID	5599/5601/5602
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around T183 & Y185 of human JNK1/2/3MAPK8 (NP_620637.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P45983/P45984/P53779
Protein Weight	35kDa/44kDa/48kDa/27kDa/52kDa
Shipping	Ice bag