

Phospho-JunB-S259 Rabbit Polyclonal Antibody For Wb - P17275

Item Number: CM0001470

Introduction

Phospho-JunB-S259 rabbit polyclonal antibody validated for western blot (WB) in human, mouse, rat. Detects transcription factor JunB phosphorylated at serine 259, ideal for AP-1 signaling pathway research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-JunB-S259 Rabbit pAb
Remarks/Alias	AP-1; Phospho-JunB-S259
Species	Human
Gene ID (Human)	3726
Gene ID	3726
Immunogen	Synthetic peptide: A phospho specific peptide corresponding to residues surrounding S259 of human Phospho-JunB-S259
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human, Mouse, Rat
SWISS	P17275
Protein Weight	36kDa
Shipping	Ice bag