

## Phospho-Lamin A/C-S22 Rabbit Pab For Wb, Ip, Elisa - P02545

Item Number: CM0001075



### Introduction

High-quality rabbit polyclonal antibody targeting phosphorylated Lamin A/C at Serine 22.

Validated for WB, IP, ELISA in human and rat samples. Ideal for nuclear lamina and chromatin organization studies.

[Learn More](#)

| Parameter        | Value                                                                                                                                  |
|------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Phospho-Lamin A/C-S22 Rabbit pAb                                                                                                       |
| Remarks/Alias    | FPL; IDC; LFP; CDDC; EMD2; FPLD; HGPS; LDP1; LMN1; LMNC; MADA; PRO1; CDCD1; CMD1A; FPLD2; LMNL1; CMT2B1; LGMD1B; Phospho-Lamin A/C-S22 |
| Species          | Human                                                                                                                                  |
| GeneID (human)   | 4000                                                                                                                                   |
| GeneID           | 4000                                                                                                                                   |
| Immunogen        | Synthetic peptide                                                                                                                      |
| Source           | Rabbit                                                                                                                                 |
| Category         | Phosphorylated Antibodies                                                                                                              |
| Application      | WB, IP, ELISA                                                                                                                          |
| Cross Reactivity | Human, Rat                                                                                                                             |
| SWISS            | P02545                                                                                                                                 |
| Protein Weight   | 74kDa                                                                                                                                  |
| Shipping         | Ice bag                                                                                                                                |