

Phospho-Lat-Y171 Rabbit Pab For Wb And Elisa - O43561

Item Number: CM0001956

Introduction

Phospho-specific rabbit polyclonal antibody targeting phosphorylated tyrosine 171 of human LAT, a key adapter in T-cell receptor signaling. Suitable for WB and ELISA applications in immunology research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-LAT-Y171 Rabbit pAb
Remarks/Alias	LAT1; pp36; IMD52; Phospho-LAT-Y171
Species	Human
GeneID (human)	27040
GeneID	27040
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around Y171 of human LAT (NP_001014987.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O43561
Protein Weight	28kDa
Shipping	Ice bag