


## Phospho-Lats1-S909/Lats2-S872 Rabbit Pab - O95835 / Q9Nrm7

Item Number: CM0000703



### Introduction

Phospho-specific rabbit polyclonal antibody recognizing LATS1 (S909) and LATS2 (S872) for WB and ELISA applications. Suitable for human samples; targets key Hippo pathway tumor suppressors.

[Learn More](#)

| Parameter        | Value                                                                        |
|------------------|------------------------------------------------------------------------------|
| Product Name     | phospho-LATS1-S909/LATS2-S872 Rabbit pAb                                     |
| Remarks/Alias    | LATS1; WARTS; wts; phospho-LATS1-S909/LATS2-S872                             |
| Species          | Human                                                                        |
| GeneID (human)   | 9113/26524                                                                   |
| GeneID           | 9113/26524                                                                   |
| Immunogen        | A synthetic phosphorylated peptide around S909 of human LATS1 (NP_004681.1). |
| Source           | Rabbit                                                                       |
| Category         | Phosphorylated Antibodies                                                    |
| Application      | WB, ELISA                                                                    |
| Cross Reactivity | Human                                                                        |
| SWISS            | O95835/Q9NRM7                                                                |
| Protein Weight   | 76kDa/126kDa                                                                 |
| Shipping         | Ice bag                                                                      |