


Phospho-Lats1-T1079+Lats2-T1041 Rabbit Monoclonal Antibody - O95835 / Q9Nrm7

Item Number: CM0001759



Introduction

Rabbit monoclonal antibody recognizing phosphorylated LATS1 (T1079) and LATS2 (T1041). Ideal for Hippo signaling studies via WB, IF/ICC, IHC-P, and ELISA. Validated in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	phospho-LATS1-T1079+LATS2-T1041 Rabbit mAb
Remarks/Alias	LATS1; WARTS; wts; phospho-LATS1-T1079+LATS2-T1041
Species	Human
GeneID (Human)	9113/26524
GeneID	9113/26524
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T1079 and T1041 of human LATS1/ LATS2 (NP_004681.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95835/Q9NRM7
Protein Weight	76kDa/126kDa
Shipping	Ice bag