


Phospho-Lck-Y394 Rabbit Pab - P06239

Item Number: CM0001510



Introduction

Rabbit polyclonal antibody detects phosphorylated LCK at Tyr394 in human and mouse. Suited for WB, IF/ICC, ELISA in T-cell receptor signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-LCK-Y394 Rabbit pAb
Remarks/Alias	LSK; YT16; IMD22; p56lck; pp58lck; Phospho-LCK-Y394
Species	Human
Gene ID (Human)	3932
Immunogen	A synthetic phosphorylated peptide around Y394 of human LCK (NP_005347.3).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P06239
Protein Weight	58kDa
Shipping	Ice bag