


Phospho-Lck-Y505 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P06239

Item Number: CM0000804



Introduction

Detects human Lck phosphorylated at Y505. Suitable for WB, IP, ELISA. Ideal for studying T-cell signaling, Lck regulation, and immune function research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Lck-Y505 Rabbit pAb
Remarks/Alias	LSK; YT16; IMD22; p56lck; pp58lck; Phospho-Lck-Y505
Species	Human
GeneID (Human)	3932
GeneID	3932
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y505 of human Lck (NP_001036236.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P06239
Protein Weight	58kDa
Shipping	Ice bag