


Phospho-Lcp2-Y145 Rabbit Pab For Wb, Elisa - Q13094

Item Number: CM0000732



Introduction

Phospho-specific rabbit polyclonal antibody targeting human LCP2 at Y145. Validated for WB and ELISA; recognizes human and mouse. Ideal for T-cell receptor signaling and immunology research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-LCP2-Y145 Rabbit pAb
Remarks/Alias	IMD81; SLP76; SLP-76; Phospho-LCP2-Y145
Species	Human
GeneID (Human)	3937
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around Y145 of human LCP2 (NP_005556.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13094
Protein Weight	60kDa
Shipping	Ice bag