


Phospho-Lyn-Y397 Rabbit Mab - P07948

Item Number: CM0000445



Introduction

Phospho-specific rabbit monoclonal antibody against Lyn (Tyr397) for WB, ELISA. Cross-reacts with mouse. Ideal for studying Lyn kinase activation in immune response, hematopoiesis, and cell signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Lyn-Y397 Rabbit mAb
Remarks/Alias	JTK8; p53Lyn; p56Lyn; Phospho-Lyn-Y396
Species	Human
GeneID (Human)	4067
GeneID	4067
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around Y396 of human Lyn (P07948).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P07948
Protein Weight	56kDa/59kDa
Shipping	Ice bag