

Phospho-Map2K3-S189 + Map2K6-S207 Rabbit Pab For Wb, Elisa - P46734 / P52564

Item Number: CM0000457

Introduction

Phospho-specific rabbit polyclonal antibody targeting phosphorylated S189 of MAP2K3 and S207 of MAP2K6, validated for Western blot and ELISA in human and mouse samples. Ideal for studying stress-activated MAP kinase signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MAP2K3-S189 + MAP2K6-S207 Rabbit pAb
Remarks/Alias	MEK3; MKK3; MAPKK3; PRKMK3; SAPKK2; SAPKK-2
Species	Human
GeneID (human)	5606/5608
GeneID	5606/5608
Immunogen	A synthetic phosphorylated peptide around S207 of human MAP2K3/MAP2K6 (NP_002747.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P46734/P52564
Shipping	Ice bag