


Phospho-Map2K4-S257 Rabbit Pab For Western Blot, Elisa - P45985

Item Number: CM0001736



Introduction

Rabbit polyclonal antibody against phosphorylated MAP2K4 at serine 257. Suitable for WB and ELISA; cross-reacts with human, mouse, rat. Key tool for studying stress kinase signaling, apoptosis, and immune homeostasis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MAP2K4-S257 Rabbit pAb
Remarks/Alias	JNKK; MEK4; MKK4; SEK1; SKK1; JNKK1; SERK1; MAPKK4; PRKMK4; SAPKK1; SAPKK-1; Phospho-MAP2K4-S257
Species	Human
Gene ID (Human)	6416
Gene ID	6416
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around S257 of human MAP2K4 (NP_003001.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45985
Protein Weight	44kDa
Shipping	Ice bag