


Phospho-Map2K4-T261 Rabbit Pab - P45985

Item Number: CM0000940



Introduction

Phospho-MAP2K4 (Thr261) rabbit polyclonal antibody for WB and ELISA. Detects human, mouse, rat MAP2K4. Suitable for studying stress-activated JNK and p38 signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MAP2K4-T261 Rabbit pAb
Remarks/Alias	JNKK; MEK4; MKK4; SEK1; SKK1; JNKK1; SERK1; MAPKK4; PRKMK4; SAPKK1; SAPKK-1; Phospho-MAP2K4-T261
Species	Human
GeneID (human)	6416
GeneID	6416
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45985
Protein Weight	44kDa
Shipping	Ice bag