

Phospho-Map2K6-S207/Map2K3-S189 Rabbit Pab For Wb - P46734 / P52564

Item Number: CM0000470

Introduction

Phospho-specific rabbit polyclonal antibody recognizing MAP2K6 (Ser207) and MAP2K3 (Ser189). Validated for Western blot and ELISA in human, mouse, and rat samples. Ideal for studying stress-activated MAP kinase signaling.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-MAP2K6-S207/MAP2K3-S189 Rabbit pAb
Remarks/Alias	MAP2K3/MAP2K6; Phospho-MAP2K6-S207/MAP2K3-S189
Species	Human
GeneID (Human)	5606 (MAP2K6) / 5608 (MAP2K3)
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46734 / P52564
Protein Weight	36kDa / 39kDa / 31kDa / 37kDa
Shipping	Ice bag