

Phospho-Map3K5 (Ask1/Mekk5) S967 Rabbit Polyclonal Antibody For Wb And Elisa - Q99683

Item Number: CM0001714



Introduction

Phospho-specific rabbit polyclonal antibody targeting human MAP3K5 (ASK1/MEKK5) phosphorylated at S967, validated for WB and ELISA applications. Useful for studying apoptosis, stress responses, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MAP3K5-S967 Rabbit pAb
Remarks/Alias	ASK1; MEKK5; MAPKKK5; Phospho-MAP3K5-S967
Species	Human
GeneID (Human)	4217
GeneID	4217
Immunogen	Synthetic peptide: A phospho specific peptide corresponding to residues surrounding S967 of human MAP3K5
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q99683
Protein Weight	155kDa
Shipping	Ice bag