

Phospho-Map3K5-S83 Rabbit Pab For Wb, Elisa - Q99683

Item Number: CM0001243



Introduction

High specificity Phospho-MAP3K5-S83 Rabbit pAb validated for Western blot (WB) and ELISA. Targets phosphorylated MAP3K5 (ASK1) at serine 83 in human. Essential for studying apoptosis, stress response and MAPK signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MAP3K5-S83 Rabbit pAb
Remarks/Alias	ASK1, MEKK5, MAPKKK5, Phospho-MAP3K5-S83
Species	Human
GeneID (Human)	4217
GeneID	4217
Immunogen	A synthetic phosphorylated peptide around S83 of human MAP3K5 (NP_005914.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q99683
Protein Weight	155kDa
Shipping	Ice bag