

Phospho-Map3K5-T845 Rabbit Polyclonal Antibody For Wb, Elisa - O35099

Item Number: CM0000853

Introduction

Phospho-MAP3K5-T845 rabbit polyclonal antibody for WB and ELISA. Detects mouse and human ASK1, a key kinase in MAPK stress and immune signaling pathways (MW ~155kDa, cross-reacts human). Ideal for apoptosis and innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MAP3K5-T845 Rabbit pAb
Remarks/Alias	ASK; ASK1; Mekk5; MAPKKK5; 7420452D20Rik; Phospho-MAP3K5-T845
Species	Mouse
GeneID (Human)	4217
GeneID (Mouse)	26408
Immunogen	Synthetic peptide / A synthetic phosphorylated peptide around T845 of mouse MAP3K5 (NP_032606.4)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O35099
Protein Weight	155kDa
Shipping	Ice bag