

Phospho-Mcl1-T163 Rabbit Pab For Wb, Elisa - Q07820

Item Number: CM0001711



Introduction

Phospho-specific rabbit polyclonal antibody targeting MCL1 phosphorylated at T163. Validated for WB and ELISA, reacts with human and mouse MCL1. Ideal for apoptosis and leukemia research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MCL1-T163 Rabbit pAb
Remarks/Alias	TM; EAT; MCL1L; MCL1S; Mcl-1; BCL2L3; MCL1-ES; bcl2-L-3; mcl1/EAT; Phospho-MCL1-T163
Species	Human
GeneID (human)	4170
GeneID	4170
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around T163 of human MCL1 (NP_068779.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q07820
Protein Weight	37kDa
Shipping	Ice bag