


Phospho-Mcm2-S139 Rabbit Polyclonal Antibody For Wb, Elisa - P49736

Item Number: CM0001331



Introduction

Phospho-MCM2-S139 Rabbit Polyclonal Antibody validated for WB, ELISA, reacts with human MCM2 phosphorylated at Ser139. Ideal for DNA replication and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MCM2-S139 Rabbit pAb
Remarks/Alias	BM28; CCNL1; CDCL1; cdc19; DFNA70; D3S3194; MITOTIN; Phospho-MCM2-S139
Species	Human
GeneID (Human)	4171
GeneID	4171
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S139 of human MCM2 (NP_004517.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P49736
Protein Weight	102kDa
Shipping	Ice bag