

Phospho-Mdm2-S260 Rabbit Pab For Wb, Elisa - Q00987

Item Number: CM0000480



Introduction

Phospho-MDM2-S260 Rabbit pAb detects human MDM2 phosphorylated at S260. Validated for WB and ELISA. Ideal for studying p53 ubiquitination, apoptosis, and cancer pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MDM2-S260 Rabbit pAb
Remarks/Alias	HDMX; LSKB; hdm2; ACTFS; Phospho-MDM2-S260
Species	Human
GeneID (Human)	4193
GeneID	4193
Immunogen	A synthetic phosphorylated peptide around S260 of human MDM2 (NP_001354919.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q00987
Protein Weight	55kDa
Shipping	Ice bag