

Phospho-Mef2C-S396 Rabbit Pab - Q06413

Item Number: CM0000617

Introduction

Phospho-MEF2C (S396) rabbit polyclonal antibody for WB and ELISA; detects human MEF2C at ~51kDa. Ideal for investigating MEF2C phosphorylation in transcription regulation, neurodevelopment, and myogenesis.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-MEF2C-S396 Rabbit pAb
Remarks/Alias	NEDHSIL; DEL5q14.3; CSDQLq14.3; Phospho-MEF2C-S396
Species	Human
GeneID (Human)	4208
GeneID	4208
Immunogen	A synthetic phosphorylated peptide around S396 of human MEF2C (NP_002388.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q06413
Protein Weight	51kDa
Shipping	Ice bag