

Phospho-Mek1 (Ps221) Rabbit Polyclonal Antibody For Wb, Elisa - Q02750

Item Number: CM0000474



Introduction

Rabbit polyclonal antibody specific for phosphorylated MEK1 at Ser221. Suitable for WB and ELISA, cross-reacts with human, mouse, and rat. Ideal for investigating MAPK/ERK signaling pathways and associated diseases.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MEK1-pS221 Rabbit pAb
Remarks/Alias	MEL; CFC3; MEK1; MKK1; MAPKK1; PRKMK1; Phospho-MEK1-pS221
Species	Human
Gene ID	5604
Immunogen	A synthetic phosphorylated peptide around S221 of human MEK1 (NP_002746.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02750
Protein Weight	43kDa
Shipping	Ice bag