

Phospho-Mek1-S217 Rabbit Pab For Wb - Q02750

Item Number: CM0000782



Introduction

Phospho-specific rabbit polyclonal antibody detecting MEK1 phosphorylated at Ser217. Validated for WB and ELISA with reactivity in human, mouse, and rat. Essential tool for MAPK/ERK pathway research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MEK1-S217 Rabbit pAb
Remarks/Alias	MEL; CFC3; MEK1; MKK1; MAPKK1; PRKMK1; Phospho-MEK1-S217
Species	Human
Gene ID (Human)	5604
Gene ID	5604
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S217 of human MEK1 (NP_002746.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02750
Protein Weight	43kDa
Shipping	Ice bag