

Phospho-Mek1-T286 Rabbit Pab - Q02750

Item Number: CM0001739



Introduction

Phospho-MEK1-T286 rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat MEK1 phosphorylated at T286. Ideal for MAPK/ERK pathway research and drug development applications.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MEK1-T286 Rabbit pAb
Remarks/Alias	MEL; CFC3; MEK1; MKK1; MAPKK1; PRKMK1; Phospho-MEK1-T286
Species	Human
GeneID (Human)	5604
GeneID	5604
Immunogen	A synthetic phosphorylated peptide around T286 of human MEK1 (NP_002746.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02750
Protein Weight	43kDa
Shipping	Ice bag