

## Phospho-Mek1-T292 Rabbit Monoclonal Antibody For Wb, Elisa - Q02750

Item Number: CM0001568



### Introduction

Phospho-MEK1-T292 rabbit monoclonal antibody validated for WB and ELISA. Detects human, mouse, rat MEK1 (MAP2K1). Essential for MAPK/ERK signaling, cell growth, differentiation. 43 kDa. Q02750.

[Learn More](#)

| Parameter        | Value                                                                                                 |
|------------------|-------------------------------------------------------------------------------------------------------|
| Product Name     | Phospho-MEK1-T292 Rabbit mAb                                                                          |
| Remarks/Alias    | MEL; CFC3; MEK1; MKK1; MAPKK1; PRKMK1; Phospho-MEK1-T292                                              |
| Species          | Human                                                                                                 |
| GeneID (Human)   | 5604                                                                                                  |
| GeneID           | 5604                                                                                                  |
| Immunogen        | Synthetic Peptide: A synthetic phosphorylated peptide around T292 of human Phospho-MEK1-T292 (Q02750) |
| Source           | Rabbit                                                                                                |
| Category         | Phosphorylated Antibodies                                                                             |
| Application      | WB, ELISA                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                     |
| SWISS            | Q02750                                                                                                |
| Protein Weight   | 43 kDa                                                                                                |
| Shipping         | Ice bag                                                                                               |