

Phospho-Mek1-T292 Rabbit Polyclonal Antibody For Wb And Elisa - Q02750

Item Number: CM0000907



Introduction

Phospho-MEK1-T292 Rabbit Polyclonal Antibody for Western blot and ELISA. Detects phosphorylated MEK1 (MAP2K1) at T292 in human and rat samples. Suitable for MAPK signaling pathway studies and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MEK1-T292 Rabbit pAb
Remarks/Alias	MEL; CFC3; MEK1; MKK1; MAPKK1; PRKMK1; Phospho-MEK1-T292
Species	Human
GeneID (human)	5604
GeneID	5604
Immunogen	A synthetic phosphorylated peptide around T292 of human MEK1 (NP_002746.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q02750
Protein Weight	43kDa
Shipping	Ice bag