

Phospho-Mek1/Mek2-S217/S221 Rabbit Mab For Wb, Ip, Elisa - Q02750/P36507

Item Number: CM0000575



Introduction

Phospho-MEK1/MEK2 (Ser217/Ser221) rabbit monoclonal antibody validated for WB, IP, ELISA. Reacts with human, mouse, rat. Ideal for studying MAPK/ERK signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MEK1/MEK2-S217/S221 Rabbit mAb
Remarks/Alias	MEK1/MEK2; Phospho-MEK1/MEK2-S217/S221
Species	Human
GeneID (Human)	5604/5605
GeneID	5604/5605
Immunogen	A synthetic phosphorylated peptide around S217 & S221 of human MEK1/MEK2 (NP_002746.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02750/P36507
Protein Weight	40kDa/43kDa/44kDa
Shipping	Ice bag