


Phospho-Mek6-S207 Rabbit Pab - P52564

Item Number: CM0001787



Introduction

Phospho-specific rabbit polyclonal antibody targeting human MEK6 phosphorylated at Serine 207, suitable for Western blot and immunofluorescence applications in human and rat samples. Ideal for stress and cytokine signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MEK6-S207 Rabbit pAb
Remarks/Alias	MEK6; MKK6; MAPKK6; PRKMK6; SAPKK3; SAPKK-3; Phospho-MEK6-S207
Species	Human
GeneID (Human)	5608
GeneID	5608
Immunogen	Synthetic peptide A phospho specific peptide corresponding to residues surrounding S207 of human MAP2K6
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC
Cross Reactivity	Human, Rat
SWISS	P52564
Protein Weight	37kDa
Shipping	Ice bag