

Phospho-Mkk3-S189+Mkk6-S207 Rabbit Mab - P46734 / P52564

Item Number: CM0001586

Introduction

Phospho-specific rabbit monoclonal antibody recognizing phosphorylated MKK3 (S189) and MKK6 (S207). Validated for WB, IHC-P, and ELISA in human, mouse, and rat samples. Ideal for stress kinase signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MKK3-S189+MKK6-S207 Rabbit mAb
Remarks/Alias	MAP2K3/MAP2K6; Phospho-MKK3-S189+MKK6-S207
Species	Human
Gene ID (Human)	5606/5608
Gene ID	5606/5608
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S207 of human MKK3/MKK6(P46734/P52564).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46734/P52564
Protein Weight	36kDa/39kDa/31kDa/37kDa
Shipping	Ice bag