

Phospho-MLKL-S345 Rabbit Mab For Wb, Ip, Elisa - Q9D2Y4

Item Number: CM0012504



Introduction

High-quality rabbit monoclonal antibody specific for mouse MLKL phosphorylated at serine 345. Ideal for Western blot, immunoprecipitation, and ELISA. Explore necroptosis signaling and TNF-induced programmed cell death.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MLKL-S345 Rabbit mAb
Remarks/Alias	9130019I15Rik
Species	Mouse
GeneID (Human)	197259
GeneID	74568
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S441 of mouse MLKL (NP_001297542.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Mouse
SWISS	Q9D2Y4
Protein Weight	54kDa
Shipping	Ice bag