

Phospho-MLkl-S358 Rabbit Pab - Q8Nb16

Item Number: CM0000546



Introduction

Phospho-specific rabbit polyclonal antibody targeting human MLKL phosphorylated at S358. Validated for Western blot and ELISA, cross-reacts with rat. Suitable for necroptosis and programmed cell death research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MLKL-S358 Rabbit pAb
Remarks/Alias	hMLKL; Phospho-MLKL-S358
Species	Human
GeneID	197259
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S358 of human MLKL (NP_689862.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8NB16
Protein Weight	54kDa
Shipping	Ice bag