


Phospho-MLKL-S441 Rabbit Pab For Wb And Elisa - Q9D2Y4

Item Number: CM0000594



Introduction

Phospho-MLKL-S441 Rabbit pAb detects mouse MLKL phosphorylated at serine 441, cross-reacting with human and rat. Suitable for WB and ELISA in necroptosis and cell death research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MLKL-S441 Rabbit pAb
Remarks/Alias	9130019I15Rik; Phospho-MLKL-S441
Species	Mouse
GeneID (Human)	197259
GeneID (Mouse)	74568
Immunogen	Synthetic peptide; A synthetic phosphorylated peptide around S441 of mouse MLKL (NP_001297542.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9D2Y4
Protein Weight	54kDa
Shipping	Ice bag