


Phospho-Mre11 (Ser676) Rabbit Pab For Wb - P49959

Item Number: CM0001997



Introduction

Phospho-MRE11 (Ser676) rabbit polyclonal antibody specifically detects human MRE11 phosphorylated at Ser676 in WB. Essential for double-strand break repair, telomere maintenance, and meiosis studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MRE11 (Ser676) Rabbit pAb
Remarks/Alias	ATLD; HNGS1; MRE11A; MRE11B; Phospho-MRE11 (Ser676)
Species	Human
GeneID (Human)	4361
GeneID	4361
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S676 of human MRE11 (NP_005582.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human
SWISS	P49959
Protein Weight	81kDa
Shipping	Ice bag