

Phospho-Mtor-S2448 Rabbit Mab For Wb, Elisa - P42345

Item Number: CM0001272



Introduction

Phospho-mTOR (Ser2448) rabbit mAb for WB, ELISA. Detects human and mouse mTOR phosphorylated at S2448, critical for mTORC1 signaling. Supports cellular metabolism and cancer research studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-mTOR-S2448 Rabbit mAb
Remarks/Alias	SKS; FRAP; FRAP1; FRAP2; RAFT1; RAPT1; Phospho-mTOR-S2448
Species	Human
Gene ID (Human)	2475
Gene ID	2475
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S2448 of human mTOR (P42345).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P42345
Protein Weight	289kDa
Shipping	Ice bag