


Phospho-MYH9-S1943 Rabbit Pab - P35579

Item Number: CM0001874



Introduction

Anti-phospho-MYH9 (S1943) rabbit polyclonal antibody validated for WB and ELISA; reacts with human, mouse, and rat myosin-9. Ideal for studying cell motility, cytokinesis, and cytoskeletal reorganization.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MYH9-S1943 Rabbit pAb
Remarks/Alias	MHA; FTNS; EPSTS; BDPLT6; DFNA17; MATINS; NMMHCA; NMHC-II-A; NMMHC-IIA; Phospho-MYH9-S1943
Species	Human
GeneID (Human)	4627
GeneID	4627
Immunogen	Synthetic peptide; A synthetic phosphorylated peptide around S1943 of human MYH9 (NP_002464.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35579
Protein Weight	227kDa
Shipping	Ice bag