

Phospho-Myl9 (Ser19) Rabbit Polyclonal Antibody For Wb And Elisa - P24844

Item Number: CM0001358

Introduction

Phospho-MYL9 (Ser19) rabbit polyclonal antibody validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat samples. Ideal for studying smooth muscle and non-muscle contractile regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MYL9-S19 Rabbit pAb
Remarks/Alias	LC20; MLC2; MRLC1; MYRL2; MLC-2C; MMIHS4; Phospho-MYL9-S19
Species	Human
Gene ID (Human)	10398
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S19 of human MYL9 (NP_006088.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24844
Protein Weight	20kDa
Shipping	Ice bag