

Phospho-Myl9-T18/S19 Rabbit Pab For Wb, Ihc-P, Elisa - P24844

Item Number: CM0000583

Introduction

Phospho-specific rabbit polyclonal antibody targeting MYL9 (Thr18/Ser19) for WB, IHC-P, ELISA. Validated in human, mouse, rat. For smooth muscle contractility and cell motility research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MYL9-T18/S19 Rabbit pAb
Remarks/Alias	LC20; MLC2; MRLC1; MYRL2; MLC-2C; MMIHS4; Phospho-MYL9-T18/S19
Species	Human
GeneID (Human)	10398
GeneID	10398
Immunogen	A synthetic phosphorylated peptide around T18 & S19 of human MYL9 (NP_006088.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24844
Protein Weight	20kDa
Shipping	Ice bag