

Phospho-Myosin Light Chain 2-S19 Rabbit Mab For Wb, Elisa - P24844

Item Number: CM0001147

Introduction

Phospho-Myosin Light Chain 2-S19 Rabbit monoclonal antibody (P24844) for WB, ELISA. Specifically detects phosphorylated S19 of human Myosin Light Chain 2 in human, mouse, rat samples. Ideal for contractility and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Myosin Light Chain 2-S19 Rabbit mAb
Remarks/Alias	LC20; MLC2; MRLC1; MYRL2; MLC-2C; MMIHS4; Phospho-Myosin Light Chain 2-S19
Species	Human
Gene ID	10398
Immunogen	A synthetic phosphorylated peptide around S19 of human Myosin Light Chain 2 (NP_006088.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24844
Protein Weight	20kDa
Shipping	Ice bag