

Phospho-Myosin Light Chain 2-T18 Rabbit Mab - P24844

Item Number: CM0001915

Introduction

Phospho-Myosin Light Chain 2 (T18) Rabbit Monoclonal Antibody for WB and ELISA. Specific for human phosphorylated MYL9/MLC2, critical for smooth muscle and nonmuscle contractility research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Myosin Light Chain 2-T18 Rabbit mAb
Remarks/Alias	LC20; MLC2; MRLC1; MYRL2; MLC-2C; MMIHS4; Phospho-Myosin Light Chain 2-T18
Species	Human
GeneID	10398
Immunogen	A synthetic phosphorylated peptide around T18 of human Myosin Light Chain 2 (NP_006088.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P24844
Protein Weight	20kDa
Shipping	Ice bag