

Phospho-Myosin Light Chain 2-T18/S19 Rabbit Mab - P24844

Item Number: CM0001200

Introduction

Phospho-specific rabbit monoclonal antibody targeting Myosin Light Chain 2 at Thr18/Ser19. Validated for WB and ELISA in human samples. Ideal for research on smooth muscle contractility, cytokinesis, and cell locomotion.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Myosin Light Chain 2-T18/S19 Rabbit mAb
Remarks/Alias	LC20; MLC2; MRLC1; MYRL2; MLC-2C; MMIHS4; Phospho-Myosin Light Chain 2-T18/S19
Species	Human
GeneID (Human)	10398
GeneID	10398
Immunogen	A synthetic phosphorylated peptide around T18/S19 of human Myosin Light Chain 2 (NP_006088.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P24844
Protein Weight	20kDa
Shipping	Ice bag