


Phospho-Na⁺/K⁺-ATPase-S16 Rabbit Pab For Wb - P05023

Item Number: CM0001352



Introduction

Rabbit polyclonal antibody detecting Na⁺/K⁺-ATPase alpha-1 subunit phosphorylated at Ser16. Suitable for WB and ELISA, cross-reacts with human, mouse, and rat. For ion transport and neurological disorder research.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-Na ⁺ /K ⁺ -ATPase-S16 Rabbit pAb
Remarks/Alias	CMT2DD; HOMGSMR2; Phospho-Na ⁺ /K ⁺ -ATPase-S16
Species	Human
Gene ID	476 (ATP1A1)
Immunogen	A synthetic phosphorylated peptide around S16 of human Na ⁺ /K ⁺ -ATPase (NP_000692.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05023
Protein Weight	113kDa
Shipping	Ice bag