

Phospho-Nf-Kb P105-S932 Rabbit Mab For Wb, Elisa - P19838

Item Number: CM0001889



Introduction

Phospho-specific rabbit monoclonal antibody targeting NF-κB p105 (Ser932) for Western blot and ELISA. Detects human, mouse, and rat NF-κB p105/p50, a key transcription factor in immune and inflammatory signaling. Ideal for pathway analysis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NF-κB p105-S932 Rabbit mAb
Remarks/Alias	KBF1; EBP-1; NF-κB; CVID12; NF-κB1; NFKB-p50; NFKappaB; NF-kappaB; NFKB-p105; NF-kappa-B1; NF-kappabeta; Phospho-NF-κB p105-S932
Species	Human
Gene ID (Human)	4790
Gene ID	4790
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S932 of human Phospho-NF-κB p105-S932 (NP_001158884.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19838
Protein Weight	105kDa
Shipping	Ice bag