


Phospho-Nf-Kb P65/Rela-S529 Rabbit Pab For Wb, Elisa - Q04206

Item Number: CM0000636



Introduction

Phospho-specific NF-kB p65 (RelA) rabbit polyclonal antibody detecting Ser529 phosphorylation; validated for WB and ELISA in human, mouse, and rat samples; for inflammation and cell signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NF-kB p65/RelA-S529 Rabbit pAb
Remarks/Alias	p65; CMCU; NFKB3; AIF3BL3; Phospho-NF-kB p65/RelA-S529
Species	Human
Gene ID (Human)	5970
Gene ID	5970
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around S529 of human NF-kB p65/RelA (NP_068810.3).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04206
Protein Weight	60kDa
Shipping	Ice bag